

CERTIFICATE OF CONFORMITY

We hereby certify that the below listed product released for delivery has been manufactured and inspected according to our approved documents and procedures. During inspection the product has fulfilled the set quality criteria and performance specifications.

Product Code: 9400100 Product Name: Finntip 200 Extended

Lot number:
19115C0

Our Quality System has been found to conform to the Quality System Standard SFS-EN ISO 9001.

Joensuu, 9th February 2022



Heidi Salokangas
Quality Manager
Laboratory Plastics Essentials
Thermo Fisher Scientific Oy