



## CERTIFICATE OF CONFORMITY


We hereby certify that the below listed product released for delivery has been manufactured and inspected according to our approved documents and procedures. During inspection the product has fulfilled the set quality criteria and performance specifications.

Product Code: 9400100      Product Name: Finntip 200 Extended

Lot number:  
21061C0

Our Quality System has been found to conform to the Quality System Standard SFS-EN ISO 9001.

Joensuu,      13th      April      2021

  
Heidi Salokangas      *Kaisa Rannasmas*  
Quality Manager  
Laboratory Plastics Essentials  
Thermo Fisher Scientific Oy