

Certificate of Conformity Thermo Scientific™ Finntip™

Product Code: 9400100
Product Name: Finntip 200 Extended
Lot Number: 22071C0

Thermo Fisher Scientific guarantees that Thermo Scientific™ Finntip™ pipette tips meet the industry's highest standards for molding precision and clarity. Our operation system is certified and conforms to ISO 9001 and ISO 14001. Finntip pipette tips are manufactured under controlled environment.

In Joensuu, 7th February 2023



Elina Maaranen
Interim Quality Manager
Thermo Fisher Scientific Oy