

Certificate of Conformity Thermo Scientific™ Finntip™

Product Code: 9400130
Product Name: Finntip 200 Extended
Lot Number: 22343C0

Thermo Fisher Scientific guarantees that Thermo Scientific™ Finntip™ pipette tips meet the industry's highest standards for molding precision and clarity. Our operation system is certified and conforms to ISO 9001 and ISO 14001. Finntip pipette tips are manufactured under controlled environment.

In Joensuu, 28th August 2023

For



Elina Maaranen
Quality Manager
Thermo Fisher Scientific Oy