

Certificate of Conformity

Finntip™ Tip Sterile

Product Code: 9400133
Product Name: Finntip 200 Extended, sterile
Lot Number: 21209C1
Exp: 2026-07

We hereby confirm that above listed high quality Finntip™ product has been manufactured and inspected according to our approved documents and procedures. During inspection the product has fulfilled the set quality criteria and performance specifications.

The sterility (SAL 10⁻⁶) is ensured by irradiating this product according to ISO 11137 guidelines.

In Joensuu, 9th February 2022



Heidi Salokangas, Quality Manager
Laboratory Plastics Essentials
Thermo Fisher Scientific Oy