

Certificate of Conformity


Finntip™ Tip Sterile

Product Code: 9400133
Product Name: Finntip 200 Extended, sterile
Lot Number: 22133C0
Exp: 2027-05

We hereby confirm that above listed high quality Finntip™ product has been manufactured and inspected according to our approved documents and procedures. During inspection the product has fulfilled the set quality criteria and performance specifications.

The sterility (SAL 10⁻⁶) is ensured by irradiating this product according to ISO 11137 guidelines.

In Joensuu, 18th May 2022



Heidi Salokangas, Quality Manager
Laboratory Plastics Essentials
Thermo Fisher Scientific Oy