


## Certificate of Conformity Thermo Scientific™ Finntip™

Product Code: 9400327  
Product Name: Finntip 10  
Lot Number: 23265A0

Thermo Fisher Scientific guarantees that Thermo Scientific™ Finntip™ pipette tips meet the industry's highest standards for molding precision and clarity. Our operation system is certified and conforms to ISO 9001 and ISO 14001. Finntip pipette tips are manufactured under controlled environment.

In Joensuu, 5th October 2023

For

  
\_\_\_\_\_  
Elina Maaranen  
Quality Manager  
Thermo Fisher Scientific Oy