



## CERTIFICATE OF CONFORMITY

We hereby certify that the below listed product released for delivery has been manufactured and inspected according to our approved documents and procedures. During inspection the product has fulfilled the set quality criteria and performance specifications.

Product Code: 9401110  
Product Name: Finntip 1000

Lot number:  
21045U0

Our Quality System has been found to conform to the Quality System Standard SFS-EN ISO 9001.

Joensuu, 27th April 2021

A handwritten signature in blue ink, appearing to read "Heidi Salokangas", written over a horizontal line.

Heidi Salokangas  
Quality Manager  
Laboratory Plastics Essentials  
Thermo Fisher Scientific Oy