



CERTIFICATE OF CONFORMITY


We hereby certify that the below listed product released for delivery has been manufactured and inspected according to our approved documents and procedures. During inspection the product has fulfilled the set quality criteria and performance specifications.

Product Code: 9401110
Product Name: Finntip 1000

Lot number:
21070U0

Our Quality System has been found to conform to the Quality System Standard SFS-EN ISO 9001.

Joensuu, 8th April 2021


Heidi Salokangas *Kaisa Rannamäki*
Quality Manager
Laboratory Plastics Essentials
Thermo Fisher Scientific Oy