

CERTIFICATE OF CONFORMITY

We hereby certify that the below listed product released for delivery has been manufactured and inspected according to our approved documents and procedures. During inspection the product has fulfilled the set quality criteria and performance specifications.

Product Code: 9401110 Product Name: Finntip 1000

Lot number:
21347F0

Our Quality System has been found to conform to the Quality System Standard SFS-EN ISO 9001.

Joensuu, 24th March 2022



Heidi Salokangas
Quality Manager
Laboratory Plastics Essentials
Thermo Fisher Scientific Oy