

Certificate of Conformity Thermo Scientific™ Finntip™

Product Code: 9401110
Product Name: Finntip 1000
Lot Number: 23251H0

Thermo Fisher Scientific guarantees that Thermo Scientific™ Finntip™ pipette tips meet the industry's highest standards for molding precision and clarity. Our operation system is certified and conforms to ISO 9001 and ISO 14001. Finntip pipette tips are manufactured under controlled environment.

In Joensuu, 30st January 2024

For



Nathan Jenkins
Quality Manager
Thermo Fisher Scientific Oy