



CERTIFICATE OF CONFORMITY


We hereby certify that the below listed product released for delivery has been manufactured and inspected according to our approved documents and procedures. During inspection the product has fulfilled the set quality criteria and performance specifications.

Product Code: 9401255
Product Name: Finntip 300

Lot number:
21025D0

Our Quality System has been found to conform to the Quality System Standard SFS-EN ISO 9001.

Joensuu, 11th March 2021


Heidi Salokangas *Kaisa Rannasmaa*
Quality Manager
Laboratory Plastics Essentials
Thermo Fisher Scientific Oy