

Certificate of Conformity Thermo Scientific™ Finntip™

Product Code: 9401255
Product Name: Finntip 300
Lot Number: 23091D0

Thermo Fisher Scientific guarantees that Thermo Scientific™ Finntip™ pipette tips meet the industry's highest standards for molding precision and clarity. Our operation system is certified and conforms to ISO 9001 and ISO 14001. Finntip pipette tips are manufactured under controlled environment.

In Joensuu, 11th July 2023

For



Elina Maaranen
Interim Quality Manager
Thermo Fisher Scientific Oy