

## Certificate of Conformity Thermo Scientific™ Finntip™

Product Code: 9401420  
Product Name: Finntip 1000 Extended  
Lot Number: 22108R0

Thermo Fisher Scientific guarantees that Thermo Scientific™ Finntip™ pipette tips meet the industry's highest standards for molding precision and clarity. Our operation system is certified and conforms to ISO 9001 and ISO 14001. Finntip pipette tips are manufactured under controlled environment.

In Joensuu, 16th January 2023



---

Elina Maaranen  
Interim Quality Manager  
Thermo Fisher Scientific Oy