

CERTIFICATE OF CONFORMITY

We hereby certify that the below listed product released for delivery has been manufactured and inspected according to our approved documents and procedures. During inspection the product has fulfilled the set quality criteria and performance specifications.

Product Code: 9402010 Product Name: Finntip 5 ml

Lot number:
22119L0

Our Quality System has been found to conform to the Quality System Standard SFS-EN ISO 9001.

Joensuu, 21st June 2022



Heidi Salokangas
Quality Manager
Laboratory Plastics Essentials
Thermo Fisher Scientific Oy