

## Certificate of Conformity

### Finntip™ Tip Sterile

Product Code: 9402073  
Product Name: Finntip 5 ml, sterile  
Lot Number: 21327M0  
Exp: 2026-11

We hereby confirm that above listed high quality Finntip™ product has been manufactured and inspected according to our approved documents and procedures. During inspection the product has fulfilled the set quality criteria and performance specifications.

The sterility (SAL 10<sup>-6</sup>) is ensured by irradiating this product according to ISO 11137 guidelines.

In Joensuu, 25th February 2022



---

Heidi Salokangas, Quality Manager  
Laboratory Plastics Essentials  
Thermo Fisher Scientific Oy