

Certificate of Conformity

Thermo Scientific™ Finntip™

Product Code: 9402171
Product Name: Finntip 10 ml
Lot Number: 24092A0

Thermo Fisher Scientific guarantees that Thermo Scientific™ Finntip™ pipette tips meet the industry's highest standards for molding precision and clarity. Our operation system is certified and conforms to ISO 9001 and ISO 14001. Finntip pipette tips are manufactured under controlled environment.

In Joensuu, 24th April 2024

For



Nathan Jenkins
Quality Manager
Thermo Fisher Scientific Oy