

## Certificate of Conformity Thermo Scientific™ Finntip™

Product Code: 9402171  
Product Name: Finntip 10 ml  
Lot Number: 24161A0

Thermo Fisher Scientific guarantees that Thermo Scientific™ Finntip™ pipette tips meet the industry's highest standards for molding precision and clarity. Our operation system is certified and conforms to ISO 9001 and ISO 14001. Finntip pipette tips are manufactured under controlled environment.

In Joensuu, 26th June 2024

For



---

Nathan Jenkins  
Quality Manager  
Thermo Fisher Scientific Oy