

**CERTIFICATE OF CONFORMITY**

We hereby certify that the below listed product released for delivery has been manufactured and inspected according to our approved documents and procedures. During inspection the product has fulfilled the set quality criteria and performance specifications.

Product Code: 9404180      Product Name: Finntip Stepper 1.25 ml

Lot number: 2117314

Our Quality System has been found to conform to the Quality System Standard SFS-EN ISO 9001.

Joensuu,      19th      April      2022



---

Heidi Salokangas  
Quality Manager  
Laboratory Plastics Essentials  
Thermo Fisher Scientific Oy