

CERTIFICATE OF CONFORMITY

FOR FINNTIP™ PIPETTE TIP

We hereby certify that the below listed product released for delivery has been manufactured and inspected according to our approved documents and procedures. During inspection the product has fulfilled the set quality criteria and performance specifications. The product marked as sterile has been irradiated by

a minimum dose of 25 kGy.

Product Code: 9404193
Product Name: Finntip Stepper 2.5 ml, sterile
Lot Number: 1808516
Expiry Date: 2023-03

Our Quality System has been found to conform to the Quality System Standard SFS-EN ISO 9001.

In Joensuu, 12th May 2022



Heidi Salokangas
Quality Manager
Laboratory Plastics Essentials
Thermo Fisher Scientific Oy