

CERTIFICATE OF CONFORMITY

FOR FINNTIP™ PIPETTE TIP

We hereby certify that the below listed product released for delivery has been manufactured and inspected according to our approved documents and procedures. During inspection the product has fulfilled the set quality criteria and performance specifications. The product marked as sterile has been irradiated by

a minimum dose of 25 kGy.

Product Code: 9404203
Product Name: Finntip Stepper 5.0 ml, sterile
Lot Number: 2300215
Expiry Date: 2028-01

Our Quality System has been found to conform to the Quality System Standard SFS-EN ISO 9001.

In Joensuu, 18th July 2023

For



Elna Maaranen
Interim Quality Manager
Laboratory Plastics Essentials
Thermo Fisher Scientific Oy