

CERTIFICATE OF CONFORMITY


FOR FINNTIP™ PIPETTE TIP

We hereby certify that the below listed product released for delivery has been manufactured and inspected according to our approved documents and procedures. During inspection the product has fulfilled the set quality criteria and performance specifications. The product marked as sterile has been irradiated by
a minimum dose of 25 kGy.

Product Code: 9404213
Product Name: Finntip Stepper 12.5 ml, sterile
Lot Number: 1812417
Expiry Date: 2023-05

Our Quality System has been found to conform to the Quality System Standard SFS-EN ISO 9001.

In Joensuu, 25th May 2021


Heidi Salokangas
Quality Manager *Kaisa Rannasmaa*
Laboratory Plastics Essentials
Thermo Fisher Scientific Oy