

CERTIFICATE OF CONFORMITY


We hereby certify that the below listed product released for delivery has been manufactured and inspected according to our approved documents and procedures. During inspection the product has fulfilled the set quality criteria and performance specifications.

Product Code: 9405050 Product Name: Finntip 1000 Wide

Lot number: 21322A2

Our Quality System has been found to conform to the Quality System Standard SFS-EN ISO 9001.

Joensuu, 20th May 2022



Heidi Salokangas
Quality Manager
Laboratory Plastics Essentials
Thermo Fisher Scientific Oy