


## Certificate of Conformity Thermo Scientific™ Finntip™

Product Code: 94060117  
Product Name: Finntip Flex 10, refill kit  
Lot Number: 22238C0

Thermo Fisher Scientific guarantees that Thermo Scientific™ Finntip™ pipette tips meet the industry's highest standards for molding precision and clarity. Our operation system is certified and conforms to ISO 9001 and ISO 14001. Finntip pipette tips are manufactured under controlled environment.

In Joensuu, 3rd October 2022

  
\_\_\_\_\_  
Elin Maaranen  
Interim Quality Manager  
Thermo Fisher Scientific Oy