



CERTIFICATE OF CONFORMITY


We hereby certify that the below listed product released for delivery has been manufactured and inspected according to our approved documents and procedures. During inspection the product has fulfilled the set quality criteria and performance specifications.

Product Code: 94060120
Product Name: Finntip Flex 10

Lot number:
20141C0

Our Quality System has been found to conform to the Quality System Standard SFS-EN ISO 9001.

Joensuu, 4th March 2021


Heidi Salokangas *Kaisa Rannamäki*
Quality Manager
Laboratory Plastics Essentials
Thermo Fisher Scientific Oy