



CERTIFICATE OF CONFORMITY

We hereby certify that the below listed product released for delivery has been manufactured and inspected according to our approved documents and procedures. During inspection the product has fulfilled the set quality criteria and performance specifications.

Product Code: 94060310
Product Name: Finntip Flex 200

Lot number:
19289A0

Our Quality System has been found to conform to the Quality System Standard SFS-EN ISO 9001.

Joensuu, 14th June 2021

A handwritten signature in blue ink, appearing to read "Heidi Salokangas".

Heidi Salokangas
Quality Manager
Laboratory Plastics Essentials
Thermo Fisher Scientific Oy