

Certificate of Conformity Thermo Scientific™ Finntip™

Product Code: 94060310
Product Name: Finntip Flex 200
Lot Number: 22228A0

Thermo Fisher Scientific guarantees that Thermo Scientific™ Finntip™ pipette tips meet the industry's highest standards for molding precision and clarity. Our operation system is certified and conforms to ISO 9001 and ISO 14001. Finntip pipette tips are manufactured under controlled environment.

In Joensuu, 14th November 2022



Elinä Maaranen
Interim Quality Manager
Thermo Fisher Scientific Oy