

CERTIFICATE OF CONFORMITY

We hereby certify that the below listed product released for delivery has been manufactured and inspected according to our approved documents and procedures. During inspection the product has fulfilled the set quality criteria and performance specifications.

Product Code: Product Name:
94060320 Finntip Flex 200

Lot number:
20165B0

Our Quality System has been found to conform to the Quality System Standard SFS-EN ISO 9001.

Joensuu, 11th August 2022



Heidi Salokangas
Quality Manager
Laboratory Plastics Essentials
Thermo Fisher Scientific Oy