

**CERTIFICATE OF CONFORMITY**


We hereby certify that the below listed product released for delivery has been manufactured and inspected according to our approved documents and procedures. During inspection the product has fulfilled the set quality criteria and performance specifications.

Product Code:    Product Name:  
94060320        Finntip Flex 200

Lot number:  
22059B0

Our Quality System has been found to conform to the Quality System Standard SFS-EN ISO 9001.

Joensuu,        13th        June        2022

  
Heidi Salokangas  
Quality Manager  
Laboratory Plastics Essentials  
Thermo Fisher Scientific Oy