

CERTIFICATE OF CONFORMITY

We hereby certify that the below listed product released for delivery has been manufactured and inspected according to our approved documents and procedures. During inspection the product has fulfilled the set quality criteria and performance specifications.

Product Code: 94060320
Product Name: Finntip Flex 200

Lot number:
22148B0

Our Quality System has been found to conform to the Quality System Standard SFS-EN ISO 9001.

Joensuu, 10th August 2022



Heidi Salokangas

Quality Manager

Laboratory Plastics Essentials

Thermo Fisher Scientific Oy