

**CERTIFICATE OF CONFORMITY**

We hereby certify that the below listed product released for delivery has been manufactured and inspected according to our approved documents and procedures. During inspection the product has fulfilled the set quality criteria and performance specifications.

Product Code: 94060510  
Product Name: Finntip Flex 300

Lot number:  
21179A0

Our Quality System has been found to conform to the Quality System Standard SFS-EN ISO 9001.

Joensuu, 4th March 2022



---

Heidi Salokangas  
Quality Manager  
Laboratory Plastics Essentials  
Thermo Fisher Scientific Oy