

**CERTIFICATE OF CONFORMITY**

We hereby certify that the below listed product released for delivery has been manufactured and inspected according to our approved documents and procedures. During inspection the product has fulfilled the set quality criteria and performance specifications.

Product Code: 94060517      Product Name: Finntip Flex 300, refill kit

Lot number:  
22195A0

Our Quality System has been found to conform to the Quality System Standard SFS-EN ISO 9001.

Joensuu,      23rd      August      2022



---

Elina Maaranen  
Interim Quality Manager  
Laboratory Plastics Essentials  
Thermo Fisher Scientific Oy