

Certificate of Conformity Thermo Scientific™ Finntip™

Product Code: 94060517
Product Name: Finntip Flex 300
Lot Number: 24070A0

Thermo Fisher Scientific guarantees that Thermo Scientific™ Finntip™ pipette tips meet the industry's highest standards for molding precision and clarity. Our operation system is certified and conforms to ISO 9001 and ISO 14001. Finntip pipette tips are manufactured under controlled environment.

In Joensuu, 22nd March 2024

For



Nathan Jenkins
Quality Manager

Thermo Fisher Scientific Oy