

CERTIFICATE OF CONFORMITY

We hereby certify that the below listed product released for delivery has been manufactured and inspected according to our approved documents and procedures. During inspection the product has fulfilled the set quality criteria and performance specifications.

Product Code: 94060520 Product Name: Finntip Flex 300 1000 / bag

Lot number:
21259A0

Our Quality System has been found to conform to the Quality System Standard SFS-EN ISO 9001.

Joensuu, 5th November 2021



Heidi Salokangas
Quality Manager
Laboratory Plastics Essentials
Thermo Fisher Scientific Oy