

CERTIFICATE OF CONFORMITY


We hereby certify that the below listed product released for delivery has been manufactured and inspected according to our approved documents and procedures. During inspection the product has fulfilled the set quality criteria and performance specifications.

Product Code: 94300110 Product Name: Labtip Yellow

Lot number:
22071D0

Our Quality System has been found to conform to the Quality System Standard SFS-EN ISO 9001.

Joensuu, 15th February 2023



Elina Maaranen
Interim Quality Manager
Laboratory Plastics Essentials
Thermo Fisher Scientific Oy