

Certificate of Quality

Thermo Scientific LabTip™


Product Code: 94300120
Product Name: Labtip Yellow
Lot Number: 24053D0
Expiry date: 2029-02

We hereby certify that the below listed product released for delivery has been manufactured and inspected according to our approved documents and procedures. During inspection the product has fulfilled the set quality criteria and performance specifications. Our operation system is certified and conforms to ISO 9001 and ISO 14001. LabTip pipette tips are manufactured under controlled environment and are certified by validation of the manufacturing process to be free of RNase, DNase and endotoxin contamination.

RNase	< 10 ⁻⁹ Kunitz units/μl
DNase	< 10 ⁻⁷ Kunitz units/μl
Endotoxins	< 0.06 EU/ml

In Joensuu, 11th July 2024

For



Nathan Jenkins
Quality Manager
Thermo Fisher Scientific Oy