

Certificate of Quality  
Thermo Scientific ClipTip™

Product Code: 94410050  
Product Name: ClipTip 384 12.5  
Lot Number: 22108B0  
Expiry date: 2027-04

Thermo Fisher Scientific guarantees that Thermo Scientific ClipTip™ Pipette tips meet the industry's highest standards for molding precision and clarity. Our operation system is certified and conforms to ISO 9001 and ISO 14001. ClipTip pipette tips are manufactured under controlled environment and are certified by validation of the manufacturing process to be free of RNase, DNase, ATP and endotoxin contamination.

|            |                                   |
|------------|-----------------------------------|
| DNA        | < 30 pg                           |
| RNase      | < $10^{-9}$ Kunitz units/ $\mu$ l |
| DNase      | < $10^{-7}$ Kunitz units/ $\mu$ l |
| Endotoxins | < 0.06 EU/ml                      |
| ATP        | < $10^{-13}$ mg/ $\mu$ l          |

In Joensuu, 27th

April 2022



---

Heidi Salokangas  
Quality Manager  
Thermo Fisher Scientific Oy