

## Certificate of Quality

### Thermo Scientific ClipTip™

Product Code: 94410060  
Product Name: ClipTip 12.5 Ext  
Lot Number: 23354C0  
Expiry date: 2028-12


Thermo Fisher Scientific guarantees that Thermo Scientific ClipTip™ Pipette tips meet the industry's highest standards for molding precision and clarity. Our operation system is certified and conforms to ISO 9001 and ISO 14001. ClipTip pipette tips are manufactured under controlled environment and are certified by validation of the manufacturing process to be free of RNase, DNase, ATP and endotoxin contamination.

DNA	< 30 pg
RNase	< 10 <sup>-9</sup> Kunitz units/μl
DNase	< 10 <sup>-7</sup> Kunitz units/μl
Endotoxins	< 0.06 EU/ml
ATP	< 10 <sup>-13</sup> mg/μl
PCR inhibitors	Not present

In Joensuu, 3rd

January 2024

For

  
Elina Maaranen

Quality Manager

Thermo Fisher Scientific Oy