


Certificate of Quality Thermo Scientific ClipTip™

Product Code: 94410070
Product Name: ClipTip 384 12.5 Ext
Lot Number: 23036C0
Expiry date: 2028-02

Thermo Fisher Scientific guarantees that Thermo Scientific ClipTip™ Pipette tips meet the industry's highest standards for molding precision and clarity. Our operation system is certified and conforms to ISO 9001 and ISO 14001. ClipTip pipette tips are manufactured under controlled environment and are certified by validation of the manufacturing process to be free of RNase, DNase, ATP and endotoxin contamination.

DNA	< 30 pg
RNase	< 10 ⁻⁹ Kunitz units/μl
DNase	< 10 ⁻⁷ Kunitz units/μl
Endotoxins	< 0.06 EU/ml
ATP	< 10 ⁻¹³ mg/μl

In Joensuu, 15th February 2023



Eina Maaranen
Interim Quality Manager
Thermo Fisher Scientific Oy