

Certificate of Quality

Thermo Scientific ClipTip™

Product Code: 94410100
Product Name: ClipTip 384 30
Lot Number: 23280D0
Expiry date: 2028-10

Thermo Fisher Scientific guarantees that Thermo Scientific ClipTip™ Pipette tips meet the industry's highest standards for molding precision and clarity. Our operation system is certified and conforms to ISO 9001 and ISO 14001. ClipTip pipette tips are manufactured under controlled environment and are certified by validation of the manufacturing process to be free of RNase, DNase, ATP and endotoxin contamination.

DNA	< 30 pg
RNase	< 10 ⁻⁹ Kunitz units/μl
DNase	< 10 ⁻⁷ Kunitz units/μl
Endotoxins	< 0.06 EU/ml
ATP	< 10 ⁻¹³ mg/μl

In Joensuu, 19th October 2023

For 

Elina Maaranen
Quality Manager
Thermo Fisher Scientific Oy