

Certificate of Quality

Thermo Scientific ClipTip™

Product Code: 94410150
Product Name: ClipTip 384, 125
Lot Number: 24053B0
Expiry date: 2029-02

Thermo Fisher Scientific guarantees that Thermo Scientific ClipTip™ Pipette tips meet the industry's highest standards for molding precision and clarity. Our operation system is certified and conforms to ISO 9001 and ISO 14001. ClipTip pipette tips are manufactured under controlled environment and are certified by validation of the manufacturing process to be free of RNase, DNase, ATP and endotoxin contamination.

DNA	< 30 pg
RNase	< 10 ⁻⁹ Kunitz units/μl
DNase	< 10 ⁻⁷ Kunitz units/μl
Endotoxins	< 0.06 EU/ml
ATP	< 10 ⁻¹³ mg/μl
PCR inhibitors	Not present

In Joensuu, 1st March 2024

For



Nathan Jenkins
Quality Manager
Thermo Fisher Scientific Oy