

Certificate of Quality
Thermo Scientific ClipTip™

Product Code: 94410217
Product Name: ClipTip 20 Reload, Low Retention
Lot Number: 22072C0
Expiry date: 2027-03

Thermo Fisher Scientific guarantees that Thermo Scientific ClipTip™ Pipette tips meet the industry's highest standards for molding precision and clarity. Our operation system is certified and conforms to ISO 9001 and ISO 14001. ClipTip pipette tips are manufactured under controlled environment and are certified by validation of the manufacturing process to be free of RNase, DNase, ATP and endotoxin contamination.

DNA	< 30 pg
RNase	< 10 ⁻⁹ Kunitz units/μl
DNase	< 10 ⁻⁷ Kunitz units/μl
Endotoxins	< 0.06 EU/ml
ATP	< 10 ⁻¹³ mg/μl

In Joensuu, 30th March 2022



Heidi Salokangas
Quality Manager
Thermo Fisher Scientific Oy