



Certificate of Quality Thermo Scientific ClipTip™

Product Code: 94410250
Product Name: ClipTip 50 , Low Retention
Lot Number: 20265C0
Expiry date: 2025-09

Thermo Fisher Scientific guarantees that Thermo Scientific ClipTip™ Pipette tips meet the industry's highest standards for molding precision and clarity. Our operation system is certified and conforms to ISO 9001, ISO 13485 and ISO 14001. ClipTip pipette tips are manufactured under controlled environment and are certified by validation of the manufacturing process to be free of RNase, DNase, ATP and endotoxin contamination.

DNA	< 30 pg
RNase	< 10 ⁻⁹ Kunitz units/μl
DNase	< 10 ⁻⁷ Kunitz units/μl
Endotoxins	< 0.06 EU/ml
ATP	< 10 ⁻¹³ mg/μl

* A detailed certificate doc. 15283950 is supplied together with each product package.

In Joensuu, 8th January 2021

Heidi Salokangas
Quality Manager
Thermo Fisher Scientific Oy