

Certificate of Quality Thermo Scientific ClipTip™

Product Code: 94410517
Product Name: ClipTip 300 Reload, Low Retention
Lot Number: 22284B0
Expiry date: 2027-10

Thermo Fisher Scientific guarantees that Thermo Scientific ClipTip™ Pipette tips meet the industry's highest standards for molding precision and clarity. Our operation system is certified and conforms to ISO 9001 and ISO 14001. ClipTip pipette tips are manufactured under controlled environment and are certified by validation of the manufacturing process to be free of RNase, DNase, ATP and endotoxin contamination.

DNA	< 30 pg
RNase	< 10 ⁻⁹ Kunitz units/μl
DNase	< 10 ⁻⁷ Kunitz units/μl
Endotoxins	< 0.06 EU/ml
ATP	< 10 ⁻¹³ mg/μl

In Joensuu, 13th October 2022



Elina Maaranen
Interim Quality Manager
Thermo Fisher Scientific Oy