



## Certificate of Quality Thermo Scientific ClipTip™

Product Code: 94410717  
Product Name: ClipTip 1000 Reload, Low Retention  
Lot Number: 21365B0  
Expiry date: 2026-12

Thermo Fisher Scientific guarantees that Thermo Scientific ClipTip™ Pipette tips meet the industry's highest standards for molding precision and clarity. Our operation system is certified and conforms to ISO 9001, ISO 13485 and ISO 14001. ClipTip pipette tips are manufactured under controlled environment and are certified by validation of the manufacturing process to be free of RNase, DNase, ATP and endotoxin contamination.

|            |                                   |
|------------|-----------------------------------|
| DNA        | < 30 pg                           |
| RNase      | < $10^{-9}$ Kunitz units/ $\mu$ l |
| DNase      | < $10^{-7}$ Kunitz units/ $\mu$ l |
| Endotoxins | < 0.06 EU/ml                      |
| ATP        | < $10^{-13}$ mg/ $\mu$ l          |

In Joensuu, 19th January 2022

  
Heidi Salokangas  
Quality Manager  
Thermo Fisher Scientific Oy