

Certificate of Quality


Thermo Scientific ClipTip™

Product Code: 94410817
Product Name: ClipTip 1250 Reload, Low Retention
Lot Number: 23157C0
Expiry date: 2028-06

Thermo Fisher Scientific guarantees that Thermo Scientific ClipTip™ Pipette tips meet the industry's highest standards for molding precision and clarity. Our operation system is certified and conforms to ISO 9001 and ISO 14001. ClipTip pipette tips are manufactured under controlled environment and are certified by validation of the manufacturing process to be free of RNase, DNase, ATP and endotoxin contamination.

DNA	< 30 pg
RNase	< 10 ⁻⁹ Kunitz units/μl
DNase	< 10 ⁻⁷ Kunitz units/μl
Endotoxins	< 0.06 EU/ml
ATP	< 10 ⁻¹³ mg/μl

In Joensuu, 14th June 2023

For 

Elina Maaranen
Interim Quality Manager
Thermo Fisher Scientific Oy