

CERTIFICATE OF CONFORMITY

We hereby certify that the below listed products are manufactured and inspected according to our approved documents and procedures before they are released for delivery. During inspection the products have to fulfill the set quality criteria and performance specifications.

Product Code:	Product Name:	Lot no:
95040470	KingFisher Flex 24DW plate	24038C2

Our Quality System has been found to conform to the Quality System Standard ISO 9001.

In Joensuu, on 24th April 2024

For



Nathan Jenkins
Quality Manager
Laboratory Plastics Essentials
Thermo Fisher Scientific Oy