




CERTIFICATE OF CONFORMITY

We hereby certify that the below listed products are manufactured and inspected according to our approved documents and procedures before they are released for delivery. During inspection the products have to fulfill the set quality criteria and performance specifications.

Product Code:	Product Name:	Lot no:	Exp:
97002610B	KF24 Flex tip comb 50 pcs/box, BC	2221010	2027-07

Our Quality System has been found to conform to the Quality System Standard ISO 9001.

In Joensuu, on 17th August 2022


Elina Maaranen
Interim Quality Manager
Laboratory Plastics Essentials
Thermo Fisher Scientific Oy