


CERTIFICATE OF CONFORMITY

We hereby certify that the below listed products are manufactured and inspected according to our approved documents and procedures before they are released for delivery. During inspection the products have to fulfill the set quality criteria and performance specifications.

Product Code:	Product Name:	Lot no:	Exp:
97002610B	KF24 Flex tip comb, BC	2406910	2029-03

Our Quality System has been found to conform to the Quality System Standard ISO 9001.

In Joensuu, on 22nd March 2024

For   
Nathan Jenkins  
Quality Manager  
Laboratory Plastics Essentials  
Thermo Fisher Scientific Oy