

CERTIFICATE OF CONFORMITY

We hereby certify that the below listed products are manufactured and inspected according to our approved documents and procedures before they are released for delivery. During inspection the products have to fulfill the set quality criteria and performance specifications.

Product Code:	Product Name:	Lot no:
97002610	KF24 Flex tip comb	2419110

Our Quality System has been found to conform to the Quality System Standard ISO 9001.

In Joensuu, on 1st August 2024

For



---

Nathan Jenkins  
Quality Manager  
Laboratory Plastics Essentials  
Thermo Fisher Scientific Oy